



# Bravo Quick Reference Guide

Operator's Manual should be consulted before using this autoclave

## BASIC OPERATION

1

2

3

**Fill – Steam Distilled Water Only**  
(if required)

Manual  
E030

Automatic  
A040

Top filling

Front filling

OR

Back

Turn unit off and then on

Reset  
(if required)

4

5

6

**Selecting a cycle**  
  
**Starting a cycle**  
  
**End of cycle**

7

8

7

**Manual Drain**  
(if required)

## KEYPAD

Reset  
Hold for 3 sec.

Interrupt  
Hold for 3 sec.

**CYCLE DESCRIPTIONS**

**Porous or Hollow**  
– Textiles, Hollow, Lumened, Solid  
(Eg.: Smocks, handpieces, scopes, Phaco)

**Solid**  
– Solid, non-lumened instrument  
(Eg.: Scissors, mirrors, pliers, scalers, burr block)

**Wrapped:**  
– Instruments that are wrapped in a cloth or paper wrap or pouched in a paper/paper or paper/plastic pouch. Also used when mixed loads are processed.

**Unwrapped:**  
– Instruments that are placed without any type of packaging.

## VENTILATION

**Front**

**Side**

## FILLING – Steam Distilled Water Only

**Manual Fill (if required)**

**Automatic Fill**

## DRAINING

**Manual Drain (if required)**

**Automatic Drain**

## MAINTENANCE SCHEDULE

To ensure trouble-free performance, operators must follow this preventative maintenance schedule.

DAILY	WEEKLY	EVERY YEAR	SPARE PARTS / ACCESSORIES
<div>Door Gasket<ul style="list-style-type: none"><li>It is recommended to check the door gasket for wear and replace if required.</li><li>Clean the door seal with a wet cloth.</li></ul></div> <div>Waste Bottle<ul style="list-style-type: none"><li>Completely empty the waste tank every time you refill the water reservoir.</li></ul></div> <div>Data Logger<ul style="list-style-type: none"><li>If using a data logger, transfer the data on the USB drive to your computer.</li></ul></div> <div>Chamber Drain Filter<ul style="list-style-type: none"><li>Inspect, clean or replace the drain filter.</li></ul></div> <div></div>	<div>Sterilization Chamber<ul style="list-style-type: none"><li>Clean the chamber and appropriate accessories with a wet cloth. If required, use a neutral detergent (a mild detergent containing no alkalies or bleaches). Carefully rinse with distilled water and remove any residue.</li><li>Clean chamber drain filter</li><li>Ensure area around Bravo is dust free</li></ul></div> <div><b>MONTHLY</b> Safety Valve<ul style="list-style-type: none"><li>Test the safety valve located on the back of the unit.</li></ul></div> <div><b>EVERY 3 – 6 MONTHS</b> Biological Filter<ul style="list-style-type: none"><li>Inspect and replace the filter, if required.</li></ul><div></div></div>	<div>Door Seal<ul style="list-style-type: none"><li>Replace every 1000 cycles or 12 months (whichever is first).</li></ul></div> <div><b>GENERAL</b> Voltage<ul style="list-style-type: none"><li>Prior to installation, ensure input voltage is greater than 220V</li></ul></div> <div>Printer Paper<ul style="list-style-type: none"><li>If printer paper is not printing properly, remove paper tighten the paper around the roll and retry.</li></ul></div>	<div><ul style="list-style-type: none"><li>Door Gasket 48000050000</li><li>Biological Filter 47200010000</li><li>Chamber Drain Filter 97290161</li><li>Thermal Paper STXX0250000</li><li>Chamber Rack 17/17V C1BP583000Y</li><li>Chamber Rack 21V C1BG534000Y</li><li>Tray, 17/17V C1XP076000Y</li><li>Tray 21V C1BG345000Y</li><li>Tray Extractor STXX0080000</li><li>Pouch Rack C1BP553000Y</li><li>Water Filling Jug A0XP0010000</li><li>Water In-Take Pull Kit 01-111774S</li><li>Direct-to-Drian Kit 01-111775S</li><li>Data Logger (G) 01-111729</li><li>Data Logger (T) 01-111730</li></ul></div>

Best Practices: SciCan recommends the final user carefully choose the most appropriate sterilization cycle according to the recommendations of their leading infection control authorities and local regulatory guidelines / recommendations.

Your Infection Control Specialist™

www.scican.com

Bravo, Dri-Tec® logo and Your Infection Control Specialist are trademarks of SciCan Ltd.

SD-346, Bravo QRG-EN